# Rheodyne® Titan**EX**®

Low-Pressure Fluidic Valves

Rheodyne TitanEX valves are a high performance, low-pressure fluidic platform with zero dead-volume. Available in low-cost, small footprint design, the TitanEX is the perfect choice when designing low-pressure instruments for the analytical, biotech, diagnostic, and industrial markets.

### Low Pressure, High Performance

The revolutionary TitanEX was designed from the ground up for low pressure (<125 psi, 9 bar) high performance fluid selection and switching applications.

### Principle of Operation

The Rheodyne TitanEX uses shear valve technology; the fluid inlet (stator) remains stationary as a grooved seal on the rotor surface rotates to change port connections. These valves are available in multiple position/port configurations.

### **Unique Tubing Connection System**

The TitanEX tubing connection system is unique, patented, and fittingless... all the parts are internal to the valve. The design reduces the tubing/valve interface to a hand-tightening operation that requires no tools. Simply insert the tubing into the ports and tighten the spanner nut. This connection system will accommodate either 1/8" or 1/16" tubing by just exchanging the internal ferrules.

# Advanced Composites Mean No Maintenance

Long-life is a principal design criteria of the TitanEX platform. The valve body and internal wear components make extensive use of modern material science by employing advanced composite polymers throughout. These highly inert and wear-resistant materials allow the valve to be actuated over the full operating temperature range without maintenance\* during its lifetime.

\*Within established ranges.



#### Low Cost, Small Footprint

With cost being a frequent OEM concern, the TitanEX is a low-cost, reliable, and quality solution in a space efficient size that can be mounted virtually anywhere.

#### Available With Integrated Driver Board

An optional driver board integrates with the TitanEX valves to provide the added functionality of motor drive and valve control without significantly affecting the remarkably small footprint of this multi-position valve. OEM customers only need to provide the digital control signals and 24V DC power in order to achieve random access actuation and position feedback. All valves may be controlled by BCD, I²C, UART, Pulse or Dual Pulse standards; two-position valves may also be controlled with level logic. In the case where multiple devices need to be controlled, I²C communication allows up to 128 devices to be connected to a single instrument. The default configuration is BCD control.

Eliminating customer development of board and firmware means shorter product development cycles and reduced time-to-market.


A Certified OEM Business Partner of ACCEPTAGE ACCEPTAGE Systec, Rheodyne, Upchurch, Isolation Technologies

info@biotech.se - www.biotech-sweden.com




## Titan**EX**™


### Dimensions given in inches and [millimeters]




6-Position, 7-Port TitanEX



6-Position, 7-Port TitanEX with PCB



10-Postion, 11-Port TitanEX with PCB



Titan	$-\mathbf{v}$	D	
			Tab II

Valves with PCB	Valves without PCB	Description	Tubing OD	Flow Passage Diameter
MLP777-601	MLP777-201	2-Position, 6-Port Injector	1/16"	0.016"
MLP777-603	MLP777-203	2-Position, 6-Port Double 3-Way	1/16"	0.016"
MLP777-605	MLP777-205	6-Position, 7-Port Selector	1/16"	0.040"
MLP777-606	MLP777-206	6-Position, 7-Port Selector	1/8"	0.040"
MLP777-612	MLP777-212	2-Position, 6-Port Injector	1/8"	0.060"
MLP777-616	MLP777-216	6-Position, 7-Port Selector	1/8"	0.060"
MLP777-624	MLP777-224	6-Position, 6-Port Switching Valve	1/8"	0.060"
MLP777-805		6-Position, 7-Port Selector, RPC-9*	1/16"	0.040"
MLP778-605	MLP778-205	10-Position, 11-Port Selector	1/16"	0.040"
MLP778-606	MLP778-206	10-Position, 11-Port Selector	1/8"	0.040"
MLP778-607	MLP778-207	2-Position, 10-Port Switching Valve	1/16"	0.040"
*Alternate mate				

TitanEX", UltraLife", and Rapid Replacement Pod" are trademarks of IDEX Health & Science LLC. DuraLife", RheBuild", and Rheodyne" are registered trademarks of IDEX Health & Science LLC. PEEK" polymer is a trademark of Victrex plc. © 2009 IDEX Health & Science LLC

R4-A4E:500-02/2009



INNOVATIVE LABORATORY PRODUCTS

A Certified OEM Business Partner of HEALTH & SCIENCE Systec, Rheodyne, Upchurch, Isolation Technologies

info@biotech.se - www.biotech-sweden.com

