



EVERYTHING FOR THE LAB

*Biotech AB – your single provider of
HPLC and fluidic scientific
customized parts & assemblies*

DEGASi[®] High Flow Degasser



*“An excellent choice for
high-throughput applications...”*



DEGASi HIGH FLOW Degasser enables efficient degassing of many types of liquids at high flow rates, up to 500ml/min.




“An excellent choice for high-throughput applications...”


APPLICATIONS

- Clinical chemistry analysers
- Dispensing applications
- Dissolution testing
- High-throughput applications
- Ink degassing

DO Removal Rate (Water, 16 kPa, 25 °C)



Pressure Drop (Water, 16 kPa, 25 °C)



GENERAL SPECIFICATIONS

Maximum Flow Rate:

500 mL/min (Water)

Maximum Channel Pressure:

0.1 MPa

Wetted Materials:

Silicone, PVC, Nitrile Rubber

Liquid Connection:

Barb for ID 8 mm

DO Removal Rate:

> 60% (300 mL/min)

Applicable Solvents:

Deionized water, pure water, neutral buffer

Dimensions:

W 95 L 315 H 140 mm

POWER REQUIREMENT

24V Power supply Included

Required input to power supply:

100-240V 50/60Hz

Part Number	Number of Channels	Internal Volume
HF.500-A	1	60 ml
HF.500-S	1	60 ml

HF.500-A is an OEM version with an open-frame structure and uses different power delivery. Please contact us for more information.



View from our office windows

BIOTECH

Biotech Sweden is a company highly specialized in chromatography accessories and consumables. Biotech is located on the Swedish Westcoast – south of Gothenburg – at the beautiful Onsala peninsula, covering 15.000 sq.m. land and 700 sq.m. production and office area.

A SIMPLE PHILOSOPHY...

We at Biotech believe that value is not "the lowest price" – our customers are the reason we exist – our customers always come first. We have not lost sight where we come from and therefore we try harder to be your partner, we're not satisfied unless you are satisfied.

Thank you for your business. I know we'll continue to grow and prosper together.

Affiliated companies: **Biotech USA LLC**
281 County Road C2W,
Roseville MN 55113, USA
bmathieu@biotechusa.us

BioNik Inc.
3397-19 Obuchi
Fuji, Shizuoka 417-0801
Japan
TEL +81-545-38-9125
info@bionikinc.com
www.bionikinc.com

We accept Visa® and MasterCard®
- no minimum order.



Biotech AB. P.O. Box 133, SE-439 23 Onsala, Sweden.
Telephone +46 (0)300-56 91 80. Fax +46 (0)300-56 91 81
info@biotech.se – www.biotech.se

